Blueprint[®] cremix

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Dansk

Blueprint[®] cremix

Dust-Free Alginate Impression Material

DeTREY Blueprint[®] cremix is a dust-free alginate impression material complying with ISO 1563 and EN 21563, respectively. The powder has been especially formulated to eliminate dust during mixing.

COMPOSITION

Diatomaceous earth (amorphous SiO₂)
Potassium alginate
Calciumsulfate dihydrate
Polypropylene glycol
Disodium orthophosphate
Potassium hexafluorotitanate
Sodium aluminosilicate (zeolith type A)
Magnesium oxide
Ultramarine blue pigment
Peppermint oil

INDICATIONS

Dental impressions for the fabrication of

- study models
- models of the opposing dentition
- working models for manufacturing of orthodontic appliances, removable prosthesis and temporary restorations.

CONTRAINDICATIONS

Do not use in patients with a known allergy to peppermint oil or other components of the material.

WARNINGS

Blueprint cremix may cause allergic reactions in patients sensitised to peppermint oil.

In isolated cases allergic reactions have been reported after impression taken with Blueprint cremix.

INTERACTIONS WITH DENTAL MATERIALS

Use of Tray Adhesive

For optimum adhesion, a tray adhesive for alginate materials, e.g. DeTREY[®] Fix, may be used.

STEP-BY-STEP INSTRUCTIONS

Proportioning

Alginate impression materials are fine powders which tend to settle during transit and storage. After transferring the contents of the bag to the plastic container, always loosen the powder before use by inverting the container twice. Sufficient volume for this has been allowed for in the container. For correct proportioning, scoop powder lightly and level off; do not tap the scoop.

Sectional impression:	1 level scoop of powder	+ lower mark of water measure
Full impression:	2 level scoops of powder	+ middle mark of water measure
Large impression:	3 level scoops of powder	+ upper mark of water measure

If a thinner or thicker mix is required use slightly more or less water.

The standard powder : water ratio¹ of 48 ml alginate powder (2 scoops) and 34 ml water (middle mark of water measure) results in 39 ml impression material.

¹ For testing according to ISO 1563: Powder : water ratio = 14.7 g powder : 34 ml water.

Mixing

Machine mixing

Mix for 8 to 12 seconds (or according to manufacturer's instructions). Machine mixing will slightly shorten both working and setting times.

Handmixing

Mix vigorously in a suitable bowl for the time stated in the table below.

Working and Setting Times²

Working and setting times can be modified by altering the temperature of the mixing water (see Table).

Water temperature (°C)	Mixing time (seconds)	Minimum working time (minutes)	Approx. setting time at 35 °C from start of mix (minutes)
15	30	2'15"	3'00"
23	30	1'30"	2'10"
31	30	1'00"	1'30"

The above setting times assume that the material is placed in the mouth 30 seconds from the end of the mix. In practice the time of placement may vary altering working and setting times. In general, at normal room temperature, the impression should remain in the mouth for approximately 1 minute.

Impression Taking

Have the patient rinse the mouth with warm water. Load and seat the tray into the mouth immediately after completion of mixing and retain, under slight pressure for approximately 1 minute. Remove the impression from the mouth and rinse free from saliva and debris, under cold running water.

Cleaning, Disinfection, Sterilisation

For disinfection immerse the impression in a 2% aqueous solution of glutaraldehyde for 10 minutes (e. g. Impresept³). After disinfection, rinse impression under running water for approximately 15 seconds.

Casting

It is recommended that the model be cast from the impression immediately after removal from the mouth. If this is not possible then wrap the impression in a damp napkin and keep in a sealed polyethylene bag to minimise any possible dimensional change.

Blueprint cremix impressions will retain their accuracy for a considerable time if stored in a humidifier maintaining at least 90% relative humidity but it is preferable that the model be cast within 48 hours of taking the impression.

Compatibility with cast materials has been confirmed for Moldano blau (Heraeus Kulzer), Vel Mix Stone (Kerr), Fuji Rock EP (GC), Calspar Rapid DP, Calstone Hard DP and Calstone Super 4.

STORAGE

Store Blueprint cremix in a tightly closed container at temperatures between 10 and 25 °C for longer periods of time.

BATCH NUMBER AND EXPIRY DATE

The batch number should be quoted in all correspondence which requires identification of the product. Do not use after expiry date.

If you have any questions, please contact:

Manufacturer:

DENTSPLY DeTREY GmbH

De-Trey-Str. 1

78467 Konstanz

GERMANY

Phone +49 (0) 75 31 5 83-0

Distributor:

DENTSPLY Limited

Building 1

Aviator Park, Station Road

Addlestone, KT15 2PG

Phone: +44 (0)19 32 85 34 22

© DENTSPLY DeTREY 2001-06-25

² The quoted working time is based on the use of distilled water but there is no significant change if tap water of any hardness is used.

³ Impresept is a trademark of 3M ESPE, Seefeld, Germany.

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DENTSPLY DeTrey GmbH